

PRO₂XY® - 5 bar

- Purity maxi. > 95%
- Oxygen available upon request or continuously
- Economic oxygen production
- Simple and reliable production technology
- Constant flow and pressure of oxygen no matter the consumption (MIL'S patent)
- Pressure at 5 bar for single stage network pressure
- PROCOM 2 control device (see page 32)
- Very compact system
- Low maintenance requirements
- Complies with ISO 10083, European Pharmacopoeia Oxygen 93 and NF S95 175




Characteristics

- Flow from 4 to 29 Nm³/h @ 93%
- Oxygen exhaust pressure : 5 bar
- O₂ measurement by heated and regulated paramagnetic sensor and with back pressure regulator
- CO / CO₂ infrared monitoring (optional)
- Dew point monitoring of O₂ (optional)
- Ambient O₂ monitoring (optional)
- Power consumption monitoring (optional)
- Working temperature from 5 °C to 45 °C
- Continuous running 24 hours (7/7)
- Power supply : 1-phase 230V - 50 / 60 Hz (100W)

PRO ₂ XY®	OX 3	OX 5	OX 7	OX 13	OX 26
Flow @ 95% (Nm ³ /h)	3	4.5	7	13	26
Flow @ 93% (Nm ³ /h)	4	5	9	14.3	28.6
Air inlet pressure	6.5 bar				
Oxygen exhaust pressure	between 4.5 and 5.5 bar				
Required air flow (Nm ³ /h)	40	60	93	150	300
Oxygen storage tank volume @ 5 bar	200 litres ▪		500 litres ▪		1000 litres ▪
Weight (kg)	290	440	640	900	1455
Required air plant in 11 bar HOSPITAIR® PACK S or G S	MVA 8	MVB 12	MVB 16	MVC 19	MVE 46

▪ Storage tank included : capacity of 200 or 500 litres mounted on the generator, capacity of 1000 litres delivered separately

Process sequence




Air is a mixture of different gases with a proportion of 78% nitrogen, 21% oxygen, and 1% argon and others gases. Oxygen Generator separates oxygen from compressed air through a pressure swing adsorption (PSA) process. The synthetic zeolite of the molecular sieve adsorbs nitrogen and concentrates the oxygen up to a level of 95%. Furthermore, the oxygen flow is kept constant thanks to the two cylinders filled with molecular sieve (synthetic zeolite).




Range complying with 2007/47/CE directive

Dimensions (mm)



	A	C
OX 3	1210	1890
OX 5	1240	2015
OX 7	1480	2140
OX 13	1580	2295



Configuration

Measurements	O ₂	O ₂ network pressure	H ₂ O O ₂	CO ₂	CO	H ₂ O air	Air pressure	Ambient O ₂	Power consumption
Single line	✓	✓	●	●	●	●	✓	●	●
Multi-line N° 1	✓	✓	●	●	●	●	✓	●	●
Multi-line N° 2	✓							●	●
Multi-line N° 3	✓							●	●

✓ standard equipment ● optional


PRO₂XY® - 5 bar - references

		OX 3	OX 5	OX 7	OX 13	OX 26
Generator	Single line	920423	920424	920425	921487	923854
Generator multi-line	Line N° 1	922311	922312	922313	922314	923855
	Line N° 2	922832	922833	922834	922835	923856
	Line N° 3	922948	922949	922950	922951	923857
1 Paramagnetic O ₂ sensor *		Foreseen 1 per generator		622272		

* Electrochemical oxygen analyser: consult us.

Optional accessory references

PROCOM 2 extension for generator single line or line N° 1 (when an option is selected)	363171
2 Manual CO / CO ₂ infrared analyser	724383
Automatic CO / CO ₂ infrared analyser	724385
3 Air hygrometry monitoring	822284
O ₂ hygrometry monitoring	823735
4 Wall-mounting analyser: O ₂ / ambient air	622709
Power consumption monitoring	see page 34
Portable ambient oxygen analyser	622765
Uninterruptible power supply (1 per generator)	624030



Specifications are subject to be changed without notice.